



# Allerbind

Panel 1 Inhalative and Food Allergens

Alder Birch Hazel Grass mix Rye Mugwort Ribwort pla Cat Horse Dog

Alternaria a Egg white Milk Peanut Hazelnut Carrot Wheal flour Soya bean

Alder Birch Hazel Oak Grass mix

Panel 4
Pediatric
Allergens

Food All

Hazelnut
Walnut
Almond
Milk
Egg white
Egg yolk
Casein
Potato
Celery
Carrot
Tomato
Cod
Crab
Orange
Apple
Wheat flor
Rye flour
Sesame

Birch Grass mix Cat Dog






Contact Allergy Proofing today

HealthPro Systems Pvt Ltd

te Tower, Mc Nichols Road, 2nd Lane, Chetpet, Chennai-600 031



Allergy Proofing - Your trusted partner in combating allergies

### **Allerbind**

Introducing one of the most extensive, precise, and detailed allergy

testing services that support allergy specialists in arriving at an accurate diagnosis of the allergen and a detailed quantitative assessment of the extent of the allergic condition with a high degree

Latest CLIA and ELISA Analyzer - intended for in vitro diagnosi

- ✓ Ability to perform both CLIA and ELISA estimations
- ✓ In-built multiple interpretation options, duplicate/ triplicate well options, and curve plotting
- ✓ A luminometer and optical system designed to perform assays in microtiter strips, read, and calculate the results

# advanced screen allergy panel

#### Introducing one of the most extensive, precise, and detailed allergy testing:



Inhalative and **Food Allergens** 



Inhalative Allergens



Panel 3 Food Allergens



Pediatric

# Redefining

- ✓ Serum-based allergy test
- ✓ Advanced immunoblot assay
- ✓ Multi-test capability
- ✓ Multi-allergen panel-based screening
- √ No initial sample load required
- ✓ Same day reporting
- Automated report generation
- ✓ High specificity
- ✓ Quantitative analysis
- ✓ In-built std O.C. for every test
- ✓ For analytics and reporting
- √ Reporting unit IU/ml

### **Allerbind**

## What's unique about Allerbind allergy testing?



A completely dedicated lab focusing exclusively on allergy testing to identify potentially harmful allergens

State-of-the-art in vitro serological allergy testing methodology for detailed analysis Comprehensive range of ELISA and CLIA tests designed to diagnose and provide detailed insights into the nature of the allergic condition

Detailed quantitative immunological profiling of the allergens and not just IgE level assessments

Consistent, automated, reproducible, accurate reports across time and geographical locations

Detailed, granular reporting that provides elaborate insights into the nature of the allergic condition, supporting allergy specialists in taking informed treatment decisions Complete line of allergens grouped into panels based on source - food allergens seasonal allergens, drug allergens, and so on





Ability to detect multiple allergies

Special pediatric panel of allergenic extracts to assess specific food allergies in children, since their food habits vary from those of adults

Home visit facilities for allergy testing so patients can get the testing done from the comfort of their homes

Rapid response time with results available online at the click of a button